



ATRAec23



Human pluripotent stem cells induced to differentiate into neural derivatives when exposed to ATRA or ec23[®]

CONSISTENT, POTENT STEM CELL DIFFERENTIATION

ec23

ec23[®] possess enhanced biological activity when compared with ATRA.

ATRA degrades in the presence of light and can cause variable results. ec23[®] is light stable and maintains its chemical stability. Choose ec23[®] and achieve the highest level of consistency in your stem cell differentiation research.

View the benefits of ec23[®]:

- Remains stable over prolonged periods, even when exposed to light
- No loss of biological activity
- Enhanced neurogenesis
- Decreased heterogeneity

To find out more visit:

visit www.reinnervate.com/ec23



reinnervate